

FOR OPTIMAL LYING COMFORT

Double rings on the 4 corners of the mattress, to loosely attach bedsheets to prevent "hammocking"

36 compartments with dynamic SmartCell air-filling

Combination of middle layer and low friction mattress cover to minimize shear forces and optimize immersion



**PRESSURE
ULCER
TREATMENT**

Optimal pressure redistribution

Pain relief

Extremely suitable for immobile users

Pain relief, unrest reduction, pressure ulcer prevention. These are just some examples of the immediate relief provided by the Vicair Mattress 415.

The design offers optimal pressure redistribution, reduction of shear forces, optimal immersion and stable – yet comfortable – positioning. Even with existing decubitus category I-IV (clinically proven).

If a situation requests an individual solution: custom sizes and shapes are available on request.

The Vicair Mattress 415 is low maintenance and comes with a fluid proof cover.

It can be used on any type of bedframe and bedframe position (e.g. fowler's position).

CPR can be performed on a patient without the use of a board.

The Vicair SmartCell air-filling makes the Vicair Mattress 415 an extremely lightweight and non-powered alternative to other air-technique replacement mattresses.

USED WITH

Pressure ulcers (category I-IV), rheumatism, burn victims, terminal illness, spasms, kyphosis, pain management.

TECHNICAL SPECIFICATIONS

Height	15 cm
Average weight	10 kg
Max. user weight	250 kg

MATTRESS SIZES

Standard

On request

Not available

Inch	width															
	27	28	30	32	34	36	38	40	42	44	46	48				
CM	68	70	75	80	85	90	95	100	105	110	115	120				
56	140															
72	180															
74	185															
76	190															
78	195															
80	200															
82	205															
84	210															
86	215															
88	220															