

LIGHTWEIGHT POSITIONING SUPPORT, COMFORTABLE AND EASY TO USE



More stability thanks to the shape of the integrated foam base

Air valves automatically adapt the cushion to your body shape and weight

Comfortable and rain proof

The Vicair Liberty Profile is an easy to use and comfortable wheelchair cushion for occasional use. Designed for comfort in combination with a low to medium level of support and pressure ulcer prevention. The Vicair Liberty Profile has an integrated foam base with a semi-circle cut out at the location of the sitting bones. The combination of SmartCells and foam-base supports pelvic stability and offers a high level of comfort in the IT-Coccyx zone. The added second valve allows an even easier automatic adaption to the body shape. This Vicair cushion is ready-for-use, rain proof and very lightweight.

PRODUCT PERFORMANCE

Skin protection / Anti decubitus	▲▲▲▲△
Pressure redistribution	▲▲▲▲△
Shear force reduction	▲▲▲▲△
Microclimate control	▲▲▲▲△
Impact & vibration damping	▲▲▲▲△

Positioning

Lateral stability	▲▲▲▲△
Frontal stability	▲▲▲▲△
Sliding reduction	▲▲▲▲△
Allows manual adjustments	▲▲▲▲△
Allows easy transfers	▲▲▲▲△

Maintenance & Safety

Durability	▲▲▲▲△
Reliability	▲▲▲▲△
Functionality when misused	▲▲▲▲▲
Low maintenance level	▲▲▲▲▲
Level of hygiene	▲▲▲▲▲

Functionality & Application

Comfort	▲▲▲▲▲
Low cushion weight	▲▲▲▲▲
Available size range	▲▲▲▲△
Maximum user weight	▲▲▲▲▲
Low seat to floor height	▲▲▲▲▲
Auto adjustments	▲▲▲▲△

USED WITH

Osteoarthritis/arthritis, rheumatism, spondylitis, pain relief, mobility problems, in a geriatric environment

TECHNICAL SPECIFICATIONS

Height	6 cm
Average weight	600 grams
Max. user weight	150 kg

CUSHION SIZES

- Standard
- Not available

		width											
		14	15	16	17	18	19	20	21	22	23	24	
		CM	35	38	40	43	45	48	50	53	55	58	60
depth	14	35											
	15	38											
	16	40											
	17	43											
	18	45											
	19	48											
	20	50											
	21	53											
	22	55											
	23	58											
	24	60											